



California Medical Systems

567 San Nicolas Dr., Ste. 280

Newport Beach, CA 92656

P: 949/719-6767

F: 949/719-6771

REAL WORLD TESTING PLAN 2023

GENERAL INFORMATION

Plan Report ID Number: Calmed_RWT_2023_5

Developer Name: California Medical Systems

Product Name(s): Practice Expert

Version Number(s): 2018

Certified Health IT Product List (CHPL) ID(s): 15.04.04.1190.Prac.18.00.1.181222

Developer Real World Testing Page URL: cal-med.com/onc.html

JUSTIFICATION FOR REAL WORLD TESTING APPROACH

California Medical System's Practice Expert v. 2018 is marketed for small to medium sized practices. The Real World Testing will focus on scenarios that occur in small to medium sized practices, encompassing single and multiple provider offices. All modules and criteria that share the purpose of Application Access will be tested. Those criteria include: § 170.315(g)(7) Application Access – Patient Selection, § 170.315(g)(8) Application Access – Data Category Request, § 170.315(g)(9) Application Access – All Category Request.

STANDARDS UPDATES (INCLUDING STANDARDS VERSION ADVANCEMENT PROCESS (SVAP) AND UNITED STATES CORE DATA FOR INTEROPERABILITY (USCDI))

Standard (and version): N/A

Update certification criteria and associated product: N/A

Health IT Module CHPL ID: N/A

Method used for standard update: N/A

Date of ONC ACB notification: N/A

Date of customer notification (SVAP only): N/A

Conformance measure: N/A

USCDI updated certification criteria (and USCDI version): N/A

Summaries

Description Summaries

Criteria	Title	Description	
§170.315(g)(7)	Application Access – Patient Selection	(i)	Receive a request to uniquely identify a patient and return a token or an ID that can be used to obtain patient data through an API.
§170.315(g)(8)	Application Access – Data Category Request	(i)(A)	Respond to API calls for specific patient data using a token or an ID to return the specified data sets for the specified patient.
§170.315(g)(9)	Application Access – All Category Request	(i)(A)	Respond to API calls for specific patient data using a token or an ID and return a CCDA file of all categories for specified patient.

Justification Summaries

California Medical Systems has a FHIR based API that allows authorized third parties to access patient data. We currently have this functionality, but as of yet, no clients or third parties have requested access. In order to test this in the real world, we will pretend to be a faux third party accessing patient data for any clients that we recruit into testing this metric. We will confirm with the client and query the database to validate that the data received via the API is accurate and in the correct format. We will de-identify the CCDA data and run it through a validator to confirm it meets the requirements.

Expected Outcomes

The Real World Testing will show that Practice Expert conforms to the following certification criteria: § 170.315(g)(7) Application Access – Patient Selection, § 170.315(g)(8) Application Access – Data Category Request, § 170.315(g)(9) Application Access – All Category Request.

REAL WORLD TESTING PLAN FIVE

CRITERIA	METRIC	METHOD	JUSTIFICATION	EXPECTED OUTCOME
170.315(g)(7) Application Access – Patient Selection	Display successful return of patient id/token when a third party requests such data.	Log the total number of patient selection requests. Also, pull patient tokens/patient IDs via Postman to confirm success.	Illustrate the volume of successful patient authentication events for patient id/token requests and whether such requests are successful.	Successful screenshots and log data showing that a patient was located and the patient id/token was sent to the third party requestor.
170.315(g)(8) Application Access – Data Category Request	Display successful return of patient data when a patient id/token and a category is requested.	Log the total number of full data category requests. Also, pull a data category request via Postman to confirm success.	Illustrate the volume of successful patient category requests and whether such requests are successful.	Successful screenshots and log data showing that a patient was located using the token/patid and that the data was sent to the third party requestor for the category requested.
170.315(g)(9) §170.315(g)(9) Application Access – All Category Request	Display a successful return of all categories for a patient when a patient id/token and all categories are requested.	Log the total number of full data requests. Also, pull a full data request via Postman to confirm success.	Illustrate the volume of successful patient all data requests and whether such requests are successful.	Successful screenshots and log data showing that a patient was located using the token/patid and that the data was sent to the third party requestor for all categories.

Care Settings

California Medical Systems markets to small and medium sized clinics, solely for ambulatory care.

Relied Upon Software

California Medical System's Practice Expert v. 2018 does not require any additional software for the criteria in Real World Testing Plan Five.

Key Milestones

Client Recruitment	12/1/2022 – 2/28/2023
Scheduling/Coordination with clients on time frame to conduct Real World Testing	1/1/2023 – 3/31/2023
Conduct Real World Testing	5/1/2023 – 8/1/2023
Parse Results	8/1/2023 – 8/31/2023
Conduct 2 nd attempt for Real World Testing in the case of anomalies or issues.	9/1/2023 – 11/1/2023
Prepare Final Results	11/1/2023 – 1/30/2024

Attestation

This Real World Testing plan is complete with all required elements, including measures that address all certification criteria and care settings. All information in this plan is up to date and fully addresses the health IT developer's Real World Testing requirements.

Authorized Representative Name: Tom Magami

Authorized Representative Email: tmagami@cal-med.com

Authorized Representative Phone: 949/719-6767

Authorized Representative Signature

A handwritten signature in cursive script, appearing to read "Tom Magami", is written over a horizontal line.

Date: 10/30/2022